ENG / UA / PL


Breathalyser

Alco Control 9.0


medica+
JAPAN TECHNOLOGY

medica+

JAPAN TECHNOLOGY

ENG

medica+
JAPAN TECHNOLOGY

Instruction manual

**Professional breathalyzer
Alco Control 9.0**

Contents

Introduction	01
Humanized safe high-tech design	01-02
Parts instruction	03
Name of parts.....	04
Operation instructio	05-14
Low battery detecion.....	15
Specification	16
Packing accessories	17
Common faults and solution.....	18-19
Unit comparison table	20
Cautions.....	21
Important notes	22
Warranty terms.....	23-24
Warning.....	25

Introduction:

Thank you for purchasing our Alco control 9.0 professional alcohol tester. Our alcohol tester adopts the UK imported latest advanced high-precision fuel cell sensor technology to produce professional-grade accurate readout in seconds. Please ensure you read the manual in its entirety before operating.

Humanized safe high-tech design

- The design and important parameters such as device sampling method, numerical accuracy, and display range refers to the EN 16280 standard drafted and officially voted by the EU in the field of breath alcohol testing devices used by the general public.
- Using high accurate and sensitive UK fuel cell sensor, more reliable and with longer life time.
- The device is designed with a breath

monitoring probe module, which monitors the continuity of breathing sampling so that the test results can more accurately reflect the human blood alcohol content (BAC).

- The device is designed with three units of measure mg/L, ‰BAC, ‰BAC, which you can easily switch between

Parts Instruction

Name of parts

1. Mouthpiece

Breath test inlet

2. LCD

Display

3. Mode Key

Unit switch and records view

4. On/Off key

ON/OFF and Test

Operation Instruction

Read the instructions carefully before use.

1. How To Use


Rinse your mouth with water before using.

Avoid smoking or drinking alcohol for at least 20 minutes before using.

2. Power ON/OFF

2.1 Power ON

Install 3 AAA alkaline batteries into the device's battery compartment, press and hold the On/Off key for 2.0 seconds to turn on and it will enter to self-checking mode automatically.

2.2 Power ON

While the device is turned on (not sampling), keep press the On/Off key for 2.0 seconds till display "OFF", then device will turn off manually.


In the test result reading and records review mode, the device will turn off automatically if left idle for 30 seconds.

3. Breath Test Instructions

3.1. Intelligent self-checking

Turn on the device, will enter into the self-checking state automatically, to reset the sensor to ensure accurate reading, the following picture will be displayed.

self checking

3.2 Breath sampling

After sensor self-checking, "BLOW" flashes for 25 seconds with beep tones and the device is ready to test, please carry out the measurement within this period, otherwise, the device will power off automatically.

Ready to blow interface

3.2 Breath sampling

When sampling, keep the mouthpiece into your mouth and blow continuously for 5 seconds, only stop blowing once the beep tone stops together with status bar full displayed. Two short beep sounds will be heard and together with vibration of the motor when the sample collection is successful, test result will be displayed accordingly.

- During breath sampling, to make the test result more accurate, please blow into the mouthpiece slowly and consistently.
- Before completing the sampling test, you can't stop blowing. "Err" will display while blowing is interrupted, it need to short press ON/OFF key again to start a new test.

Err (blowing interrupted)

test result display

- At the test result or "Err" display interface, can short press ON/OFF key for a new test. And if keep pressing ON/OFF key for 2.0s, the device will power off.

3.3 Unit switch

At the test result display interface, short press the Mode key to switch the measurement units between mg/L, ‰BAC, %BAC. The next test result will be shown with the unit you selected.

‰BAC

%BAC

mg/L

Tip: If the customer requires the display of a fixed unit before leaving the factory, the unit switching function is invalid.

3.4 Test result show

If the test result exceeds the set alert level, the alarm indicator light turns on red and flashes 10 times together with "di-di" sound prompt.

While test result lower than set alarm level, the alert indicator light turns on green.

The reading range is 0.00‰BAC -4.00‰BAC.

"4.00‰ BAC" will be shown if the reading is higher than 4.00‰ BAC.

"0.00‰BAC" will be shown if the reading is lower than 0.05‰ BAC.

Take alert level setted as 0.50‰BAC for example, test result will show as following pictures:

< 0.50‰BAC

≥ 0.50‰BAC

3.5 Set alarm value

Alarm value with factory default setting (Generally there is no need to set up again). If you need to change the alarm value, in the self-checking mode, press ON/OFF key twice to enter into the alarm value setting mode, and short press ON/OFF key to select the alarm value as b2 ,b3,b5,b8, after selected short press ON/OFF key to start a new test.

Corresponding critical value and meaning are as following:

4 configurable alarm level values.

B2

$C < 0.20\text{‰BAC}$ or 0.020‰BAC or 0.100mg/L

no alarm

$C \geq 0.20\text{‰BAC}$ or 0.020‰BAC or 0.100mg/L

alarm

B3

$C < 0.30\text{‰BAC}$ or 0.030‰BAC or 0.150mg/L

no alarm

$C \geq 0.30\text{‰BAC}$ or 0.030‰BAC or 0.150mg/L

alarm

B5

$C < 0.50\text{‰BAC}$ or 0.050‰BAC or 0.250mg/L

no alarm

$C \geq 0.50\text{‰BAC}$ or 0.050‰BAC or 0.250mg/L

alarm

B8

$C < 0.80\text{‰BAC}$ or 0.080‰BAC or 0.400mg/L

no alarm


$C \geq 0.80\text{‰BAC}$ or 0.080‰BAC or 0.400mg/L

alarm

Tip: If customer has no special requirements, the factory default setting of alarm value should be used.

3.6 Test records review

Device will record the 20 most recent testing results. Hold key for 2.0s to enter into record memory mode. After entering, press the key repeatedly to review the last stored records (max. 20 records), the times of record and test result will flash on the display alternately. For example, if the 1st record is 0.50‰ BAC, then "-01-" and the BAC reading "0.50‰ BAC" flashes alternately; if the 2nd has not been tested, then the recorded data will show ".." (As per the following pictures):

Low battery detection

When there is a low battery level, "Lob" flashes with sound prompt as following picture. After "Lob" flashes 10 times the device will power off automatically, please replace with new AAA alkaline batteries.

Low battery prompt

Warning

While low battery the device will not work normally, you need to replace with new batteries, or it may not power on.

Note:

When the batteries are taken out, the stored records will be cleared.

Specifications

Sensor: UK fuel cell sensor

Response Time: 2-3s

Blowing Time: About 5s

Detection Range:

0.00 ~4.00‰BAC;

0.000~0.400%BAC;

0.000~2.000mg/L;

(If the conversion factor between ‰BAC units and mg/L units is 2.0, then 0.50 ‰BAC =

0.050 %BAC = 0.250 mg/L)

Sensitivity: 0.001% BAC

Sensor Accuracy:

+/-0.006 at 0.050% BAC;

+/-0.008 at 0.080% BAC;

>0.080% BAC: +/-10%

Breath test monitoring

Breathing interruption microphone detection

Record: 20 test records be stored

Optimal Operating Temperature: 0 °C - 50 °C

Storage Temperature: -20 °C - 70 °C

Display: About 2.0 inch LCD(blue film)

Power Supply: 3 × AAA alkaline batteries

Dimensions: 132.5 * 68.5 * 26mm

Net Weight: 87g

Packing accessories

Packing accessories may vary by different countries or operators, please check the packing list according to the order contract.

Standard packing:

- 1* Alco control 9.0 equipment
- 1* user manual
- 1* color box
- 10* Mouthpieces

Optional accessories:

- 3* AAA alkaline batteries, Mouthpieces.

Tip: If find that the device is damaged or missing accessories, please contact your supplier.

Common faults and solutions

Can't power on

1. Low battery. - Replace new batteries.
2. Batteries are upside down. - Re-install the batteries in right way.
3. Haven't held press ON/OFF key for 2.0s - Hold down the ON/OFF key for 2.0s to turn on the device.

Err

1. During breath sampling, blowing is interrupted midway. - Re-test, at breath sampling, keep blow continuously for 5s.

Two test results differ greatly in a short time

1. After drinking alcohol, the waiting time is not long enough before the alcohol test. - To get a more accurate test result, please abstain from smoking and drinking for at least 20 minutes before testing.
2. Before the test, no water or non-alcoholic

beverages was used to rinse the mouth and a large amount of alcohol residue remained in the mouth. -

To ensure the test result more closer to human Blood Alcohol Concentration, Please Rinse out your mouth with water or non-alcoholic drinks before sampling test.

Unit comparison table

NO.	Blood Alcohol Content (BAC)		Breath Alcohol Content (BrAC)
	%BAC	‰BAC	mg/L
1	0.000	0.00	0.000
2	0.010	0.10	0.050
3	0.020	0.20	0.100
4	0.030	0.30	0.150
5	0.040	0.40	0.200
6	0.050	0.50	0.250
7	0.060	0.60	0.300
8	0.070	0.70	0.350
9	0.080	0.80	0.400
10	0.090	0.90	0.450
11	0.100	1.00	0.500
12	0.110	1.10	0.550
13	0.120	1.20	0.600
14	0.130	1.30	0.650
15	0.140	1.40	0.700
16	0.150	1.50	0.750
17	0.160	1.60	0.800
18	0.170	1.70	0.850
19	0.180	1.80	0.900
20	0.190	1.90	0.950
21	0.200	2.00	1.000
...
41	0.400	4.00	2.000

Cautions

- Different people have different tolerances for alcohol which is affected by factors such as heredity, physique, gender, and weight. Therefore, different people can drink the same amount of alcohol, and the test results will be different.
- After drinking alcohol, there is a process of absorption and decomposition of alcohol in the blood, so the alcohol content in your blood can fluctuate within a few minutes.
- Store the device in its packing box, avoid exposure to sunlight and humidity.

Important notes

- Don't allow alcohol or any other liquid enter into the device, or it will cause permanent damage to the sensor.
- The frequency of use and BAC test concentration of device, will affect the sensitivity and life of the sensor. The more frequently used and higher alcohol concentration tested the sensor sensitivity will be reduced and service life will be shortened.
- Test results are for reference only, do not use this device to be 100% sure if you can drive or not. The test results may be inadmissible in court.
- For health and safety of you and others, please don't drink before drive

Warranty Terms

Dear customers,

Please keep this instruction manual properly to use the warranty card. The warranty is provided for one year from the date of purchase (or for 2 or 3 years if the extended warranty is purchased) in case of malfunction under normal operating conditions.

The warranty does not cover the following cases:

1. If the product is damaged due to improper use or lack of proper maintenance.
2. If the product is damaged as a result of disassembly or modifications.
3. If parts not supplied by our company were applied to and used with the device.
4. Damage caused by force majeure circumstances.

According to the terms, technical maintenance provided outside the warranty is subject to additional charges.

Parts of the product replaced during the warranty period remain the property of MEDICA+.


Warning

Neither the re-seller nor manufacturer accepts any liability arising from the use of the unit. Don't drink before you drive!

Note: All illustrations on this instruction are effect drawings. Manufacturer reserves the right to modify the contents of this manual without prior notice.

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